


<b>Product name:</b>		<b>StressVue: Stress Exercise Testing System</b>	
<b>Identification code(s)</b>		<b>860295 StressVue System</b>	
<b>Total weight (in Kg)</b>		<b>Total = 3.9 Kg for Wired or Wireless System</b> <b>989803142811 Direct Wire ECG Transmitter = 0.14 Kg</b> <b>989803145561 USB Connector (Direct Wired) = 0.22 Kg</b>  <b>989803145531 Wireless 2500 ECG Transmitter= 0.11 Kg</b> <b>989803145551 USB Connector (Wireless 2500)= 0.29 Kg</b>  <b>989803142811 ECG Printer = 3.51 Kg</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems	
	Address:	3000 Minuteman Road, Andover Ma	
	Zip code:	01810	
	Country:	USA	

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	Special tools required for disassembly: None	
	Distribution to EU is limited since only English is available  The StressVue system uses several different electronic assemblies that can only be used as part of the total system. The Philips branded parts that are covered by the WEEE Directive are listed above.	See photo's on pages 2 and 3

<b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
		
	None	

<b>Batteries</b>	<b>Type:</b>	<b>Location</b>
	AA Alkaline battery in Wireless Transmitter only	See page 2
<b>To be Removed</b>		

<b>Hazardous</b>	<b>Substances:</b>	<b>Location</b>
		
	None	
<b>To be Removed</b>		

Printed copies are uncontrolled

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled




**Direct Wired ECG Module with USB Connector: No battery included**



**Wireless ECG Module with USB Connector: Includes AA battery see battery door on back of ECG transmitter for removal.**

Printed copies are uncontrolled

Printed copies are uncontrolled



**Printer Used with StressVue System**