





| | | |
|-------------------------------|-------------------------|-------------------------|
| Product name: | Allura XPER FD10 | |
| Identification code(s) | 4522 901 04430 | |
| Total weight (in Kg) | 2000 Kg (approx) | |
| Producer /Manufacturer | Name company: | Philips Medical Systems |
| | Address: | Veenpluis 6 |
| | Zip code: | 5684 PC Best |
| | Country: | The Netherlands |
| | Electronic info: | Allura_FD10_579.doc |

| Recycle info | Substances: | Location |
|---|---|------------------|
| Special attention  | Items: | Location |
| | Be aware of possible contaminated system parts and materials. | See Figure 1 |
| | For dismantling activities Treatment Facilities must consider the national requirements. | |
| | For personnel that can come in contact with contaminated material, preventive measurements pursuant to national requirements must taken into account. | |
| | Removal of units/weights can cause the system to tilt! | |
| | Removal of units/weights can cause unexpected movements. | |
| | High voltage parts are marked with | |
| Hazardous  | Lead (Pb) for X-Ray Shielding | |
| | Lead (Pb) for soldering | All electronics |
| | Beryllium copper (BeCu) | Some electronics |
| | Mercury (Hg) in specific LCD Screens | LCD Screens |
| | | |
| To be removed | | |
|  Fluids / Gases | Items: | Location |
| | Transformer Oil, type | |

Note: To facilitate recycling, all plastic parts > 50 grams are marked according to ISO 11469 & ISO 1043

| | | | | | |
|---|---------------------|-----------------------------|-------------------|---|--|
| Name: WEEE Passport Allura XPER FD10 | | 12NC: 4522 901 04430 | | Version: 00 Doc.Date: 2006-09-01 | |
| Author: W.Bleyen | Pag.: 1 of 3 | Group No: 579 | Section: - | Status: Approved | |

Figure 1



For Cabinets :
 See R-Cabinet R-Passport
 See M-Cabinet R-Passport
 See X-Ray Generator Passport for Velara CV Cabinet

| | | | | | |
|--|--|--------------------------------|----------------------|---|-------------------------|
| Name: WEEE Passport Allura XPER FD10 | | 12NC: 4522 901 04430 | | Version: 00 Doc.Date: 2006-09-01 | |
| Author: W.Bleyen | | Pag.: 2 of 3 | Group No: 579 | Section: - | Status: Approved |

Recycling Passport Revision History

| Rev | Description | Date |
|-----|----------------|---------------|
| 1 | First Passport | Sept. 1, 2006 |

Recycling Passport Approval Signatures

| |
|---------------------------------------|
| Written By PMS W.Bleyen |
|---------------------------------------|

| |
|------------------------|
| Reviewed By (optional) |
|------------------------|

| |
|-------------------------------|
| Approved By PMS |
|-------------------------------|

| |
|-------------|
| Approved By |
|-------------|

| | | |
|---|-----------------------------|--|
| Name: WEEE Passport Allura XPER FD10 | 12NC: 4522 901 04430 | Version: 00 Doc.Date: 2006-09-01 |
| Author: W.Bleyen | Pag.: 3 of 3 | Group No: 579 Section: - Status: Approved |