





<b>FCC Cabinet HFO</b>	
Identification code(s)	9896 030 1360x - 9896 030 1657x - 9896 030 2069x
Total weight [kg]	± 170 kg excluding CDAS and RMMU
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	

## Recycling information

Special attention	Items	Location
	CDAS (4522 132 0893x) or CDAS ASSY II (4598 000 1447x) are not part of this passport	Figure 2
	RMMU not part of this passport	Figure 2

Batteries	Items, type	Location
	RMMU, see Magnet HFO (9896 030 1586x)	Figure 2

Hazardous substances	Items	Location
	BeCu EMC-springs on edge of Filter case	Figure 1
	BeCu EMC-springs on edge of DFU	Figure 1

Fluids and gases	Items	Location
	Not applicable	

Name:	<b>Recycling passport FCC Cabinet HFO</b>	12NC:	<b>9896 030 13600</b>	Version: <b>01</b>
Author:		Pag.: <b>1 of 2</b>	Group No: 5790	Doc.Date: <b>2010-02-12</b>
			Section: -	Status: <b>Authorized</b>

## Locations as mentioned in the Passport (Information in pictures)




Figure 1: FCC seen from Scanner room side




Figure 2: FCC seen from Technical room side

Name:	<b>Recycling passport FCC Cabinet HFO</b>	12NC:	<b>9896 030 13600</b>	Version: <b>01</b>
Author:		Group No: 5790	Section: -	Doc.Date: <b>2010-02-12</b>
	Pag.: <b>2 of 2</b>			Status: <b>Authorized</b>