




Panorama	Panorama 1.0T Covers
Identification code(s)	<p>BOX 1 COVERS HFO 4522-132-5100x consiting: FRONT COVER HFO ASSY 4522-132-3403 TOP PLATE FRONT HFO ASSY 4522-132-3428</p> <p>BOX 2 COVERS HFO 4522-132-5101x consiting: REAR COVER HFO ASSY 4522-132-3407 TOP PLATE REAR HFO ASSY 4522-132-3430</p> <p>BOX 3 COVERS HFO 4522-132-5102x consiting: SIDE COVER RGT HFO ASSY 4522-132-3416 TOPFRAME GRP ASSY 4522-132-4687 QBC C-COVER FL-RT ASSY 4522-132-5155 QBC C-COVER FT-RL ASSY 4522-132-5156</p> <p>BOX 4 COVERS HFO 4522-132-5103x consiting: SIDE COVER LFT HFO ASSY 4522-132-3414 FRONT EDGE COVER HFO 4522-132-3410 REAR EDGE COVER HFO 4522-132-34123 POST MID COVER 4522-132-3432 RG COVER LEFT 4522-132-3434 RG COVER RIGHT 4522-132-3436</p> <p>BOX 5 COVERS HFO 4522-132-5104x consiting: TOP COVER LEFT HFO ASSY 422-132-34182 TOP COVER RIGHT HFO ASSY 4522-132-3424</p> <p>COVER SCU PLEAIDES 4522 132 1701x</p>
Total weight [kg]	± 200 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	System_CoversPanorama1T_452213251000_579_v02.doc


Recycling information

Special attention	Items	Location
	Products are made of composite materials consisting of glass, fibers and resins.	
	Microphone PWB in front cover	Figure 1

Name:	Recycling Passport Panorama 1.0T Covers	12NC:	4522 132 51000	Version: 02
Author: Jos van Vliet		Pag.: 1 of 4	Group No: 5790	Doc.Date: 2007-08-01
			Section: -	Status: Authorized

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Not applicable	

Fluids and gases	Items	Location
	Not applicable	

Name: Recycling Passport Panorama 1.0T Covers		12NC: 4522 132 51000		Version: 02 Doc.Date: 2007-08-01	
Author: Jos van Vliet		Pag.: 2 of 4	Group No: 5790	Section: -	Status: Authorized

Overview of covers:



figure 1

Printed circuit board locations
(inside of covers)

Name: Recycling Passport Panorama 1.0T Covers		12NC: 4522 132 51000		Version: 02 Doc.Date: 2007-08-01	
Author: Jos van Vliet		Pag.: 3 of 4	Group No: 5790	Section: -	Status: Authorized
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007				Template version 2.	

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	8-9-2006
2	Authorised version	11-9-2006

Recycling Passport Approval Signatures

<p>Written By J. van Vliet</p>	<p>Reviewed By (optional) G.J. Plattel</p>
<p>Approved By J. van Vliet</p>	<p>Approved By A. Klaver</p>

End of document.

Name: Recycling Passport Panorama 1.0T Covers		12NC: 4522 132 51000		Version: 02 Doc.Date: 2007-08-01	
Author: Jos van Vliet		Pag.: 4 of 4	Group No: 5790	Section: -	Status: Authorized