





Printed copies are uncontrolled

<b>Product name:</b>	<b>DigiTrak+ Holter Recorder</b>	
<b>Identification code(s)</b>	<b>860277 Philips Ordering</b> <b>203005-03: part # for 24 Hour Recording Capability</b> <b>203005-04: part # for 48 Hour Recording Capability</b>	
<b>Total weight (in Kg)</b>	<b>0.2 Kg with AA Alkaline Battery</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems
	Address:	3000 Minuteman Road, Andover Ma
	Zip code:	01810
	Country:	USA

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Special Tools:</b> Need TORX 5 P wrench to remove back door	N/A
<b>Fluids / Gases</b> 	<b>Items:</b>	N/A
	None	
<b>Batteries</b> 	<b>Type:</b>	See page 3
	Alkaline 1 AA per unit	
	1 Lithium Ion button cell in unit	
<b>To be Removed</b>		
<b>Hazardous</b> 	<b>Substances:</b>	See page 2
	LCD Display contains (Hg).	
<b>To be Removed</b>		

Show locations of materials mentioned on the previous sheet. (Picture information)

Printed copies are uncontrolled



**LCD Display Location**

Printed copies are uncontrolled



**Lithium Ion  
Button Cell**



**AA Alkaline Battery Location**