





Printed copies are uncontrolled

<b>Product name:</b>	<b>8002 Compressor/Cooler Cabinet (Panorama 0.6T MRI System)</b>	
<b>Identification code(s)</b>	<b>612-909 /4510 000 02021</b>	
<b>Total weight (in Kg)</b>	<b>420</b>	
<b>Producer/ Manufacturer</b>	Name company:	Philips Medical Systems MR Technologies Finland
	Address:	Äyritie 4
	Zip code:	01511 VANTAA
	Country:	FINLAND

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	High pressure 1.65MPa helium gas inside the adsorber and the Supply and Return Flex Lines of the Compressor Unit / Use He-gas purging valve to release the pressure slowly.	Figure 3 Location 1
<b>Fluids / Gases</b> 	Oil	Figure 3 Location 2
<b>Batteries</b> 	<b>Type:</b>	<b>Location</b>
	NONE	
<b>Hazardous</b> 	<b>Substances:</b>	<b>Location</b>
	NONE	

Show locations of materials mentioned on the previous sheet. (Picture information)


**Figure 1**



**Figure 2**



**Figure 3** Adsorber inside the compressor



Location 1

Location 2

Title: Recycling Passport Panorama 0.6T Cooler Compressor  
Recycling Passport Number: 987046 Rev: A

Printed copies are uncontrolled

### Recycling Passport Revision History

Rev	Description	Date
A	Initial version	9.8.2006

### Recycling Passport Approval Signatures

Written By  Tero Virta
------------------------------

Reviewed By (optional)  Review: 989775-A
--

Approved By  Tatu Laine
-------------------------------

Printed copies are uncontrolled