




Product name	Philips Pediatric coil
45356734154	4522 131 78020 Pediatric coil 1.5T
Total weight [kg]	9 Kg
Manufacturer	Philips Medical Systems
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFCoils_Philips_PediatricCoils.doc

Recycling information

Special Attention 	Items	Location
	None	Figure

Batteries 	Items, type	Location
	None	Figure

Hazardous substances 	Items	Location
	printed circuit boards inside the coil	Figures next page
	Beryllium copper contact pins on connectors	Figures next page

Fluids and gases 	Items	Location
	None	

Name Recycling passport Philips Pediatric coil		12NC: 4522 131 78020	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben	Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



Figure 1

Name Recycling passport Philips Pediatric coil	12NC: 4522 131 78020	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 2 of 3	Group No: 5790 Section: -
		Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

Written By Paul Ebben	Reviewed By (optional) Geert-Jan Platel Harrie Schellekens
Approved By Paul Ebben	Approved By (optional)

End of document.

Name Recycling passport Philips Pediatric coil	12NC: 4522 131 78020	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -	Status: Authorized