




<b>Product name:</b>	<b>BT sensor set</b>		
<b>Identification code(s)</b>	<b>4522 132 24540</b>		
<b>Total weight (in Kg)</b>	<b>2</b>		
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Nederland B.V.	
	Address:	Veenpluis 4-6	
	Zip code:	5680 DA BEST	
	Country:	The Netherlands	
	Electronic info:	BasicTriggeringSensorSet_452213224540_579_v01.doc	


BT Sensor Set contains loose items (sensors eo)


<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Items:</b>	<b>Location</b>
	NA	

 <b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	LI ion battery	See photo next page

<b>Hazardous</b>   <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	Printed Circuit Board inside battery charger	See photo next page

Name: <b>Recycling Passport BT sensor set</b>		12NC: <b>4522 132 24540</b>		Version: <b>01</b> Doc.Date: <b>2007-05-11</b>	
Author: <b>J. Schoonderbeek</b>		Pag.: <b>1 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>

 <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	NA	



LI-ion battery



Printed Circuit Board

### Recycling Passport Revision History

Rev	Description	Date
1	First Passport	June 21, 2006

### Recycling Passport Approval Signatures

Written By

**H. Schoonderbeek**

Reviewed By (optional)

Approved By

**H. Schoonderbeek**

Approved By

Name: <b>Recycling Passport BT sensor set</b>		12NC: <b>4522 132 24540</b>	Version: <b>01</b> Doc.Date: <b>2007-05-11</b>
Author: <b>J. Schoonderbeek</b>	Pag.: <b>2 of 2</b>	Group No: <b>5790</b>	Section: <b>-</b>
			Status: <b>Authorized</b>