




Product name:	FASTTRAK HFO		
Identification code(s)	9896 030 1595 x (where x can be any digit from 1 .. 9)		
Total weight (in Kg)	285		
Producer /Manufacturer	Name company:	Philips Medical Systems Nederland B.V.	
	Address:	Veenpluis 4-6	
	Zip code:	5680 DA BEST	
	Country:	The Netherlands	
	Electronic info:	Patient_FastTrakHFO_989603015950_579_v01.doc	



Recycle info	Substances:	Location
Special attention 	Items:	Location
	Inside the triangular section of the tabletop, there is a cable and Printed Circuit Boards.	1
	This product contain a scissor mechanism and spindle. Care shall be taken during disassembly to prevent injury due to collaps of the mechanism.	2
	This product contains Printed Circuit Boards inside the indicated housings.	3
	This product contains the 'Basic Triggering Unit'. Refer to the passport of 4522 132 2456 x.	4
	All aluminium parts of this product are hexavalent chrominated.	
	This product has internal electric cabling. Disassemble covers to get access to all cabling.	

Batteries	Type:	Location
 To be Removed	This product contains lead acid batteries.	5

Hazardous		
 To be removed	The Printed Circuit Boards contain lead based solder.	1, 3

Name: Recycling Passport FastTrak Trolley		12NC: 9896 030 1595 0		Version: 00 Doc.Date: 2006-09-14	
Author: G.J. Lenting		Pag.: 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Pictures information)



Name: Recycling Passport FastTrak Trolley		12NC: 9896 030 1595 0	Version: 00 Doc.Date: 2006-09-14
Author: G.J. Lenting	Pag.: 2 of 3	Group No: 5790	Section: - Status: Authorized

Recycling Passport Revision History

Rev	Description	Date
00	First issue	September 14, 2006

Recycling Passport Approval Signatures

Written By GJ Lenting
GJ Lenting

Reviewed By H Boonen
H Boonen

Reviewed By J Bregman
J Bregman

Approved By W Evers
W Evers

Name: Recycling Passport FastTrak Trolley			12NC: 9896 030 1595 0			Version: 00 Doc.Date: 2006-09-14			
Author: G.J. Lenting		Pag.: 3 of 3		Group No: 5790		Section: -		Status: Authorized	