




<b>Product name</b>	<b>Philips Cylindrical QBC's</b>
Identification code(s)	4522 131 77010 QBC 1.5T → 35 Kg 4522 131 79010 QBC 1.5T → 15 Kg 4535 670 13080 QBC 3.0T → 85 Kg 4535 671 40260 QBC 3.0T (without booster iron) → 65 Kg
Total weight [kg]	See list above
<b>Manufacturer</b>	<b>Philips Medical Systems</b>
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFCoils_Philips_QBC.doc

### Recycling information

Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure

Hazardous substances	Items	Location
	printed circuit boards inside the coil	Figure 1
	Beryllium copper contact pins on connectors	Figure 1
	Beryllium oxide in 'load' resistor	Figure 2

Fluids and gases	Items	Location
	None	

Name <b>Recycling passport Philips QBC's</b>		12NC: <b>4522 131 77010</b>		Version: <b>00</b> Doc.Date: <b>2006-08-04</b>	
Author: <b>Paul Ebben</b>		Page <b>1</b> of <b>4</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>



Name <b>Recycling passport Philips QBC's</b>	12NC: <b>4522 131 77010</b>	Version: <b>00</b> Doc.Date: <b>2006-08-04</b>	
Author: <b>Paul Ebben</b>	Page <b>2</b> of <b>4</b>	Group No: <b>5790</b> Section: -	Status: <b>Authorized</b>

**Locations as mentioned in the Passport (Information in pictures)**

*Place sketch (es) or photo(s) here, and indicate issues using icons:*



Figure 1




Figure 2 (cover already dismantled)

Name <b>Recycling passport Philips QBC's</b>		12NC: <b>4522 131 77010</b>		Version: <b>00</b> Doc.Date: <b>2006-08-04</b>	
Author: <b>Paul Ebben</b>		Page <b>3</b> of <b>4</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>

### Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

### Recycling Passport Approval Signatures

<p><b>Written By</b> Paul Ebben</p>	<p><b>Reviewed By (optional)</b> Geert-Jan Platel Harrie Schellekens</p>
<p><b>Approved By</b> Paul Ebben</p>	<p><b>Approved By (optional)</b></p>

End of document.

Name <b>Recycling passport Philips QBC's</b>	12NC: <b>4522 131 77010</b>	Version: <b>00</b> Doc.Date: <b>2006-08-04</b>
Author: <b>Paul Ebben</b>	Page <b>4</b> of <b>4</b>	Group No: <b>5790</b> Section: -
		Status: <b>Authorized</b>