





|                        |   |
|------------------------|---|
| <b>Product name</b>    | <b>Philips Body coils</b>                           |
| Identification code(s) | 4510 000 41530 BodySp-M<br>4510 000 41540 BodySp-XL |
| Total weight [kg]      | BodySp-M → 6.2 Kg<br>BodySp-XL → 8.3 Kg             |
| <b>Manufacturer</b>    | <b>Philips Medical Systems</b>                      |
| Address                | PL 185 / Äyritie 4                                  |
| Zip code               | 01511 VANTAA  |
| Country                | FINLAND   |
| Electronic info        | RFCoils_Philips_BodyCoils.doc                       |

### Recycling information

| Special Attention   | Items | Location |
|---|-------|----------|
|  | None  | Figure   |
|   |       |          |
|   |       |          |
|   |       |          |

| Batteries   | Items, type | Location |
|---|-------------|----------|
|  | None        | Figure   |
|   |             |          |
|   |             |          |
|   |             |          |


| Hazardous substances  | Items  | Location         |
|---|--|------------------|
|  | Several printed circuit boards inside the coil | Figure next page |
|   | Beryllium copper contact pins on connectors    | Figure next page |
|   |  |                  |
|   |  |                  |

| Fluids and gases  | Items | Location |
|---|-------|----------|
|  | None  |          |
|   |       |          |
|   |       |          |
|   |       |          |

|  |  |                                |                       |   |                              |
|--|--|--------------------------------|-----------------------|---|------------------------------|
| Name<br><b>Recycling passport Philips Body coils</b> |  | 12NC:<br><b>4510 000 41530</b> |                       | Version: <b>00</b><br>Doc.Date: <b>2006-08-04</b> |                              |
| Author: <b>Paul Ebben</b>                            |  | Page <b>1</b> of <b>3</b>      | Group No: <b>5790</b> | Section: -  | Status:<br><b>Authorized</b> |

## Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



|  |                                |   |                              |
|--|--------------------------------|---|------------------------------|
| Name<br><b>Recycling passport Philips Body coils</b> | 12NC:<br><b>4510 000 41530</b> | Version: <b>00</b><br>Doc.Date: <b>2006-08-04</b> |                              |
| Author: <b>Paul Ebben</b>                            | Page <b>2</b> of <b>3</b>      | Group No: <b>5790</b> Section: -                  | Status:<br><b>Authorized</b> |

## Recycling Passport Revision History

| Revision | Description     | Date       |
|----------|-----------------|------------|
| 00       | Initial version | 2006-08-04 |

## Recycling Passport Approval Signatures

|                                  |  |
|----------------------------------|--|
| <b>Written By</b><br>Paul Ebben  | <b>Reviewed By (optional)</b><br>Geert-Jan Platel<br>Olli Friman |
| <b>Approved By</b><br>Paul Ebben | <b>Approved By (optional)</b>                                    |

End of document.

|  |                                |   |                              |
|--|--------------------------------|---|------------------------------|
| Name<br><b>Recycling passport Philips Body coils</b> | 12NC:<br><b>4510 000 41530</b> | Version: <b>00</b><br>Doc.Date: <b>2006-08-04</b> |                              |
| Author: <b>Paul Ebben</b>                            | Page <b>3</b> of <b>3</b>      | Group No: <b>5790</b> Section: -                  | Status:<br><b>Authorized</b> |