




Product name	RF Amplifier, 15 KW, 64 Mhz
Identification code(s)	4522 131 4582x
Total weight [kg]	270 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFamplifierS23_452213145820_579_v01.doc

Recycling information

Special Attention	Items	Location
	None	NA

Batteries	Items, type	Location
	Lithium Ion Energy Cell – Dallas Semiconductor 1244Y	Figure 1

Hazardous substances	Items	Location
	RF Deck Printed Circuit Boards Containing Tin / Lead Solder	Figure(s) 1,2,3,4
	Power Supply Printed Circuit Boards Containing Tin / Lead Solder	Figure 5
	Beryllium Copper Gaskets	Figure(s) 2,3

Fluids and gases	Items	Location
	Capacitors Containing Approximately 160 Milliliters Of Mineral Oil	Figure 5

Name: RF Amplifier, 15 KW, 64 Mhz		12NC: 4522 131 4582x	Version: 00
Author: Michael Loucks		Group No: 5790	Doc.Date: 2006-08-03
Page 1 of 6	Section: -	Status: Auth	

Locations as mentioned in the Passport (Information in pictures)




Figure 1
RF Deck – Front View (Power Supply Below Not Show)



Name: RF Amplifier, 15 KW, 64 Mhz		12NC: 4522 131 4582x	Version: 00 Doc.Date: 2006-08-03
Author: Michael Loucks	Page 2 of 6	Group No: 5790	Section: - Status: Auth




Figure 2
RF Deck Left Side



Note #1
Printed Circuit Board
Inside Enclosure



Note #2
Printed Circuit Board
Opposite Side



Note #3
Printed Circuit Board
Inside Tube

Name: RF Amplifier, 15 KW, 64 Mhz		12NC: 4522 131 4582x	Version: 00 Doc.Date: 2006-08-03
Author: Michael Loucks	Page 3 of 6	Group No: 5790	Section: - Status: Auth




Figure 3
RF Deck – Top View Sides Open





Figure 4
RF Deck Front Pannel



Note #4
Printed Circuit Board



Note #5
Printed Circuit Board Inside Canister




Figure 5
Power Supply Top / Front View (RF Deck Removed)



Note #6
Capacitors Located In Right Rear
Of Enclosure



Name: RF Amplifier, 15 KW, 64 Mhz		12NC: 4522 131 4582x	Version: 00
Author: Michael Loucks		Group No: 5790	Doc.Date: 2006-08-03
Page 5 of 6		Section: -	Status: Auth

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2006-07-03

Recycling Passport Approval Signatures

Written By Michael Loucks	Reviewed By (optional) Fred Saar
Approved By Rick Carpenter	Approved By Don Fowler

End of document.

Name: RF Amplifier, 15 KW, 64 Mhz		12NC: 4522 131 4582x	Version: 00
Author: Michael Loucks		Page 6 of 6	Doc.Date: 2006-08-03
		Group No: 5790	Section: -
		Status: Auth	