


Product name	Philips Flex coils
Identification code(s)	4522 132 31150 Flex-L 1.5T 4522 132 19160 Flex-L 3.0T 4522 131 81410 Flex-L 1.5T IT 4522 132 19100 Flex-M 1.5T (old) 4522 132 31170 Flex-M 1.5T 4522 132 19260 Flex-M 3.0T 4522 131 81420 Flex-M 1.5T IT 4522 132 31190 Flex-S 1.5T 4522 132 60250 Flex-S 3.0T 4522 131 81400 Flex-R 1.5T IT 4522 131 66940 C1-II 4522 131 15260 C3 4522 131 15270 C4 4522 131 66950 E1-II 4522 131 66960 R1 4522 131 15900 Cable trap 1.5T
Total weight [kg]	Flex-S → 3 KG Flex-M/R → 3.5 Kg Flex-L → 3.5 Kg C1, C3, C4 → 0,2 Kg E1, R1 → 1Kg Cable trap → 2Kg
Manufacturer	Philips Medical Systems
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFCoils_Philips_FlexCoils.doc

Recycling information

Hazardous substances	Items	Location
	Several PWB's. Some flexible some rigid.	Figures next page
	Beryllium copper contact pins on connectors	Figures next page


Name Recycling passport Philips Flex coils	12NC: 4522 132 31150	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 1 of 3	Group No: 5790 Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



Equal for all coils (coil, driver box, connector)




Flex-L coil



Flex-M coil




Flex-S coil



Flex-R IT & Flex-L IT coils



Cable trap used on C1, C3, C4, E1 and R1 coils



C1, C3 and C4 coils



E1 coil



R1 coil

Name Recycling passport Philips Flex coils		12NC: 4522 132 31150	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 2 of 3	Group No: 5790	Section: -
			Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

Written By Paul Ebben	Reviewed By (optional) Geert-Jan Platel
Approved By Paul Ebben	Approved By (optional)

End of document.

Name Recycling passport Philips Flex coils	12NC: 4522 132 31150	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -	Status: Authorized