




Product name	Invivo 4-channel Breast Coils
Identification code(s)	4522 150 40350 4-Breast 1.0T 4522 150 39190 4-Breast 1.5T 9896 030 13360 4-Breast 3.0T
Total weight [kg]	8 Kg.
Manufacturer	Invivo, corporation. USA
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_Breast-4chCoils.doc

Recycling information

Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure

Hazardous substances	Items	Location
	Several printed circuit boards inside the coil	Figure 1
	Beryllium copper contact pins on connectors	Figure 1 + 2

Fluids and gases	Items	Location
	Plastic bottle with copper sulfate	Figure 3
	0.008moles/liter CuSO ₄ and 0.034moles/liter NaCl	

Name: Recycling passport Invivo 4-channel Breast coils		12NC: 4522 150 40350		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:




Figure 1




Figure 2



Figure 3

Name: Recycling passport Invivo 4-channel Breast coils		12NC: 4522 150 40350	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 2 of 3	Group No: 5790	Section: -
			Status: Authorized	

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

Written By Paul Ebben	Reviewed By (optional) Geert Jan Platel Bert Jansen
Approved By Paul Ebben	Approved By

End of document.

Name: Recycling passport Invivo 4-channel Breast coils	12NC: 4522 150 40350	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -	Status: Authorized