



<b>Product name</b>	<b>Invivo Shoulder Coils</b>
Identification code(s)	4510 000 49230 4-shoulder 1.0T 4522 150 40200 4-shoulder 1.5T 9896 030 13400 4-shoulder 3.0T 4535-301-04121 SENSE Shoulder Coil, 1.5T 4535-301-04151 SENSE Shoulder Coil, 3.0T
Total weight [kg]	8.5 Kg.(4-ch shoulder coil) 6 Kg (8-ch shoulder coil)
<b>Manufacturer</b>	<b>Invivo, corporation. USA</b>
Address	3545 SW 47th Avenue
Zip code	Gainesville, FL 32608
Country	USA
Electronic info	RFCoils_Invivo_ShoulderCoils v01.doc

### Recycling information

<b>Special Attention</b>	<b>Items</b>	<b>Location</b>
	None	Figure

<b>Batteries</b>	<b>Items, type</b>	<b>Location</b>
	None	Figure

<b>Hazardous substances</b>	<b>Items</b>	<b>Location</b>
	Several printed circuit boards inside the coil	Figure 1,3
	Beryllium copper contact pins on connectors	Figure 1,2,4

Fluids and gases	Items	Location
	4-ch coils; Plastic bottle with copper sulfate 0.008moles/liter CuSO <sub>4</sub> and 0.034moles/liter NaCl	Figure 5
	8-ch coils; Plastic bottle with Nickel Chloride 31.5mM/L of NiCL <sub>2</sub> -6H <sub>2</sub> O	Figure 6

Locations as mentioned in the Passport (Information in pictures)



Figure 1




Figure 2




Figure 3



Figure 4




Figure 5



Figure 6