





| | |
|------------------------|--|
| Product name | Standard cabinet 33HU + 2HU drawers. |
| Identification code(s) | 4522 131 7443x |
| Total weight [kg] | 123 kg |
| Manufacturer | Philips Medical Systems Nederland B.V. |
| Address | Veenpluis 4-6 Best |
| Zip code | 5680DA |
| Country | The Netherlands |
| Electronic info | System_StandardCabinetNTDAC_452213174430_579_v02.doc |

Recycling information

| | | |
|--|---|------------------|
| Special Attention  | Items | Location |
| | Not applicable | |
| | | |
| | | |
| Batteries  | Items, type | Location |
| | Not applicable | |
| | | |
| | | |
| Hazardous substances  | Items | Location |
| | Rivet nut M12 (Zink yellow passivated) | Figure 2 D + 4 D |
| | Rivet nut M6 (Zink yellow passivated) | Figure 2 B |
| | Rivet nut M8 (Zink yellow passivated) | Figure 2 A |
| | Remark: Nuts can only be removed with e.g. a drill. | |
| | Levellers (adjustable feet) (Zink yellow passivated bolt) | Figure 2 C |
| Fluids and gases  | Items | Location |
| | Not applicable | |
| | | |
| | | |

| | | | |
|--|---------------------------|--------------------------------|---|
| Name: Recycling passport Standard cabinet 33HU + 2HU | | 12NC: 4522 131 74430 | Version: 02 Doc.Date: 2007-07-31 |
| Author: G. Geboers | Page 1 of 2 | Group No: 5790 | Section: - Status: Auth. |

Locations as mentioned in the Passport (Information in pictures)



Figure 1




Figure 2

A (4 x 4)

B (1 + 1)

C (4x)

D (4 x)



Figure 3



Figure 4

D (4 x)

Recycling Passport Revision History

| Revision | Description | Date |
|----------|-----------------|------------|
| 1 | Initial version | 2006-08-18 |

Recycling Passport Approval Signatures

| | |
|---|-------------------------------|
| Written By G. Geboers | Reviewed By (optional) |
| Approved By Geert-Jan Plattèl | Approved By |

End of document.

| | | | |
|--|---------------------------|--------------------------------|---|
| Name: Recycling passport Standard cabinet 33HU + 2HU | | 12NC: 4522 131 74430 | Version: 02 Doc.Date: 2007-07-31 |
| Author: G. Geboers | Page 2 of 2 | Group No: 5790 | Section: - |
| © PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007 | | | Status: Auth. |