





Product name	Standard cabinet 33HE + 2HE lade.
Identification code(s)	4522 131 7443x
Total weight [kg]	123 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	DasNTDACCabinet_452213174430_579_v01.doc

Recycling information

Special Attention 	Items	Location
	Not applicable	
Batteries 	Items, type	Location
	Not applicable	
Hazardous substances 	Items	Location
	Rivet nut M12 (Zink yellow passivated)	Figure 2 D + 4 D
	Rivet nut M6 (Zink yellow passivated)	Figure 2 B
	Rivet nut M8 (Zink yellow passivated)	Figure 2 A
	Remark: Nuts can only be removed with e.g. a drill.	
	Levellers (adjustable feet) (Zink yellow passivated bolt)	Figure 2 C
Fluids and gases 	Items	Location
	Not applicable	

Name: Recycling passport Stand.Cab. 33HE + 2HE Lade		12NC: 4522 131 7443 x	Version: 00 Doc.Date: 2006-08-18
Author: G. Geboers	Page 1 of 2	Group No: 5790	Section: - Status: Auth.

Locations as mentioned in the Passport (Information in pictures)



Figure 1




Figure 2

A (4 x 4)

B (1 + 1)

C (4x)

D (4 x)




Figure 3



Figure 4

D (4 x)

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2006-08-18

Recycling Passport Approval Signatures

Written By G. Geboers	Reviewed By (optional)
Approved By Geert-Jan Plattèl	Approved By

End of document.

Name: Recycling passport Stand.Cab. 33HE + 2HE Lade		12NC: 4522 131 7443 x	Version: 00 Doc.Date: 2006-08-18
Author: G. Geboers	Page 2 of 2	Group No: 5790	Section: -
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2006			Status: Auth.