




Achieva	Achieva 1.5T/3.0T Covers
Identification code(s)	<p>SET COVERS 1.5T :</p> <p>SET FRONT COVERS WE 4522 132 2755x</p> <p>SET REAR COVERS WE 4522 132 2756x</p> <p>SET TOP COVERS WE 4522 132 2757x</p> <p>SIDE COVER WE ASSY (2X) 4522 132 2678x</p> <p>SET COVERS 3.0T :</p> <p>SET FRONT COVERS 3.0T WE 4522 132 2758x</p> <p>SET REAR COVERS 3.0T WE 4522 132 2759x</p> <p>SET TOP COVERS 3.0T WE 4522 132 2780x</p> <p>SIDE COVER 3.0T WE ASSY 4522 132 2716x</p> <p>SIDE COVER WE ASSY (2X) 4522 132 2678x</p>
Total weight [kg]	± 150 kg
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	SystemCoversAchieva_452213227550_579_v01.doc


Recycling information

Special attention	Items	Location
	Products are made of composite materials consisting of glass, fibers and resins.	
	Microphone printed circuit board in front cover	Figure 1

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Not applicable	

Name: SET FRONT COVERS WE	12NC: 4522 132 27550	Version: 00 Doc.Date: 2006-09-11
Author: J. van Vliet	Page 1 of 3	Group No: 5790 Section: - Status: Auth.

Fluids and gases	Items	Location
	Not applicable	

Overview of covers:



Printed circuit board location
(inside of cover)

figure 1

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	8-9-2006
2	Authorised version	11-9-2006

Recycling Passport Approval Signatures

Written By J. van Vliet	Reviewed By (optional) G.J. Plattel
Approved By J. van Vliet	Approved By A. Klaver

End of document.

Name: SET FRONT COVERS WE	12NC: 4522 132 27550	Version: 00 Doc.Date: 2006-09-11
Author: J. van Vliet	Page 3 of 3	Group No: 5790 Section: -
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2006		Status: Auth.