




Product name	Recycling Passport LCC
Identification code(s)	98960301390x 98960301391x
Total weight [kg]	380 kg
Manufacturer	Glen Dimplex Deutschland GmbH
Address	Am Goldenen Feld 18
Zip code	95326 Kulmbach
Country	Germany
Electronic info	Gradient_LiquidCoolingCabinet_989603013900_579_v02.doc

Recycling information


Special Attention	Items	Location
	Ensure unit is drained of cooling fluids	Tank, Liquid circuits

Batteries	Items, type	Location
	Not applicable	

Fluids and gases	Items	Location
	Dowtherm SR-1 in secondary GA – circuit. Properties: Ethylene glycol 50%/50% with additives, manufacturer Dow Chemical	Tank, Fluid lines, pump Secondary
	Water with Biozid and Inhibitor in secondary GC – circuit. Properties: Biozid: 10 vol-%, Typ NX1164, supplier CWG GmbH Inhibitor: 20 vol-%, Typ AZ8104, supplier CWG GmbH	Tank, Fluid lines, pump Secondary



Name: Recycling Passport LCC		12NC: 9896 030 13900		Version: 02 Doc.Date: 2007-07-06	
Author: R. Vloet		Pag.: 1 of 4	Group No: 5790	Section: -	Status: Authorized

	Items	Location
	Residual primary cooling fluid: water with 0 to 39 vol-%(ethylene) glycol and possibly other unknown additives	Fluid lines Primary

Hazardous substances	Items	Location
	Lead is used in 3-way-valve actuator	Actuator
	PBB and PBDE are elements of PVC in edge guard	Edge guard

Name: Recycling Passport LCC		12NC: 9896 030 13900		Version: 02 Doc.Date: 2007-07-06	
Author: R. Vloet		Pag.: 2 of 4	Group No: 5790	Section: -	Status: Authorized


Locations as mentioned in the Passport



Position of Printed Circuit Board inside LCC:



Printed Circuit Board



Recycling Passport Revision History

Revision	Description	Date
1	Initial version	June 22, 2007
2	Authorized version	July 06, 2007

Recycling Passport Approval Signatures

Written By Rolf Sebastian, Development Manager GDD Robert Vloet	Reviewed By (optional) Robert Vloet Geert-Jan Plattèl
Approved By Robert Vloet	Approved By Geert-Jan Plattèl

<END OF DOCUMENT>

Name: Recycling Passport LCC	12NC: 9896 030 13900	Version: 02 Doc.Date: 2007-07-06		
Author: R. Vloet	Pag.: 4 of 4	Group No: 5790	Section: -	Status: Authorized