




Product name	Philips Head coils
Identification code(s)	4522 131 39000 Head coil 1.5T 4522 132 19830 T/R Head 3.0T 4522 131 23340 H-Head Left (PICU) II 1.5T 4522 131 23260 H-Head Right (PICU) II 1.5T
Total weight [kg]	7.5 Kg
Manufacturer	Philips Medical Systems
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFCoils_Philips_Headcoils.doc

Recycling information


Special Attention	Items	Location
	None	Figure

Batteries	Items, type	Location
	None	Figure

Hazardous substances	Items	Location
	Several printed circuit boards inside the coils.	Figures next page
	Beryllium copper contact pins on connectors	Figures next page

Name Recycling passport Philips head coils	12NC: 4522 131 39000	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 1 of 3	Group No: 5790 Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)



Head coil 1.5T



T/R Head coil 3.0T



H-Head

Name Recycling passport Philips head coils		12NC: 4522 131 39000		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 2 of 3	Group No: 5790	Section: -	Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

<p>Written By Paul Ebben</p>	<p>Reviewed By (optional) Geert-Jan Platel Harrie Schellekens</p>
<p>Approved By Paul Ebben</p>	<p>Approved By (optional)</p>

End of document.

Name Recycling passport Philips head coils	12NC: 4522 131 39000	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -
		Status: Authorized