





<b>Product name</b>	<b>Gradient Coils: TNF2, TNF3, Triton, AG30S</b>
Identification code(s)	4522 131 8002.x, 4522 131 8009.x, 4522 131 8006.x, 4522 131 8007.x
Total weight [kg]	850 kg
<b>Manufacturer</b>	<b>Tesla Engineering Limited</b>
Address	Water Lane, Storrington, Sussex
Zip code	RH20 3EA
Country	UK
Electronic info	GradientCoil_TNF2_TNF3_Triton_AG30S_452213180020_579_v01.doc

## Recycling information


<b>Special Attention</b> 	<b>Items</b>	<b>Location</b>
	Coil contains cooling water + additives, see below.	Figure 1
<b>Batteries</b> 	<b>Items, type</b>	<b>Location</b>
	No batteries	
<b>Hazardous substances</b> 	<b>Items</b>	<b>Location</b>
	Solder in electrical connections contains lead (Pb).	Figure 1,2
<b>Fluids and gases</b> 	<b>Items</b>	<b>Location</b>
	Additives to cooling water: Betz Dearborn Inhibitor AZ8104 ; 0.2mL/L Betz Dearborn Biocide Spectrus NX1100, 0.09mL/L	Figure 1

Name: <b>Recycling passport TNF2,TNF3, Triton, AG30S</b>		12NC: <b>4522 131 8002.x</b>		Version: <b>00</b> Doc.Date: <b>2006-07-27</b>	
Author: <b>J.Konijn</b>		Page <b>1</b> of <b>3</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>

## Locations as mentioned in the Passport (Information in pictures)

Place sketch(es) or photo(s) here, and indicate issues using icons:

Figure 1.




Water contains Inhibitor and Biocide



Soldered connections contain Pb (many places internal in coil)



Figure 2.



Name: <b>Recycling passport TNF2,TNF3, Triton, AG30S</b>		12NC: <b>4522 131 8002.x</b>	Version: <b>00</b> Doc.Date: <b>2006-07-27</b>
Author: <b>J.Konijn</b>	Page <b>2</b> of <b>3</b>	Group No: <b>5790</b>	Section: <b>-</b>
			Status: <b>Authorized</b>

## Recycling Passport Revision History

Revision	Description	Date
1	Initial version	27-07-2006

## Recycling Passport Approval Signatures

<b>Written By</b> Jan Konijn	<b>Reviewed By (optional)</b> G-J Plattèl
<b>Approved By</b> J.Konijn (PMS)	<b>Approved By</b> P.Smithers (Tesla)

End of document.

Name: <b>Recycling passport TNF2,TNF3, Triton, AG30S</b>	12NC: <b>4522 131 8002.x</b>	Version: <b>00</b> Doc.Date: <b>2006-07-27</b>	
Author: <b>J.Konijn</b>	Page <b>3</b> of <b>3</b>	Group No: <b>5790</b> Section: -	Status: <b>Authorized</b>