



Product name	Evic Endo Coil
Identification code(s)	4522 131 15010 Endo coil 4522 131 68330 Endo coil small
Total weight [kg]	1 Kg.
Manufacturer	Evic/Neways
Address	Voltaweg 10
Zip code	6101 XK Echt
Country	Netherlands
Electronic info	RFCoils_EvicEndocoil.doc

Recycling information

Special Attention 	Items	Location
	None	Figure

Batteries 	Items, type	Location
	None	Figure

Hazardous substances 	Items	Location
	Several printed circuit boards inside the coil.	Figure 1
	Beryllium copper contact pins on connectors	Figure 1

Fluids and gases 	Items	Location

Name: Recycling passport Evic Endo coil		12NC: 4522 131 15010		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:




Figure 1

Name: Recycling passport Evic Endo coil		12NC: 4522 131 15010		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 2 of 3	Group No: 5790	Section: -	Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

Written By Paul Ebben	Reviewed By (optional) Geert Jan Platel Harrie Schellekens
Approved By Paul Ebben	Approved By

End of document.

Name: Recycling passport Evic Endo coil	12NC: 4522 131 15010	Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -	Status: Authorized