


<b>Product name:</b>	PFEI-C-30T	
<b>Identification code(s)</b>	4522-132-2025x	
<b>Total weight (in Kg)</b>	20 Kg	
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Best
	Address:	Veenpluis 4-6
	Zip code:	5680DA Best
	Country:	The Netherlands
	Electronic info:	PFEI-C-30T

<b>Recycle info</b>	<b>Substances:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Items:</b>	
	No special attention needed.	
 <b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	No batteries present.	
<b>Hazardous</b>   <b>To be removed</b>	<b>Examples</b>	<b>Location</b>
	Depending on the implementation four or six Printed circuit board assemblies	Figure 1a and 1b
	Each Printed circuit board assembly contains 10 times 14mm beryllium copper sticky fingers.	Figure 2
 <b>Fluids / Gases</b>	<b>Items:</b>	<b>Location</b>
	No fluids and/or gases present.	


Name: <b>Recycling Passport PFEI-C-30T</b>		12NC: <b>4522 132 20250</b>		Version: <b>04</b>	
Author: <b>R. Kunnen</b>		Pag.: <b>1 of 3</b>		Doc.Date: <b>2007-08-08</b>	
		Group No: <b>5790</b>		Section: <b>-</b>	
				Status: <b>Authorized</b>	

© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007 Template version 2.


**Locations as mentioned in the Passport ( Pictures information)**



*Figure 1a, PFEI-C-30T right-side*



*Figure 1b, PFEI-C-30T left-side*



*Figure 2, Beryllium Sticky-Fingers on Printed circuit board assembly*


Name: <b>Recycling Passport PFEI-C-30T</b>		12NC: <b>4522 132 20250</b>		Version: <b>04</b> Doc.Date: <b>2007-08-08</b>	
Author: <b>R. Kunnen</b>		Pag.: <b>2 of 3</b>	Group No: <b>5790</b>	Section: <b>-</b>	Status: <b>Authorized</b>
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2007					Template version 2.

### Recycling Passport Revision History

Rev	Description	Date
1	First Passport (Draft)	March 31, 2006
2	Second Version (Concept)	27 July, 2006
3	Third Version (Authorised)	28 July, 2006
4	Fourth Version (Authorised)	08 August 2007


### Recycling Passport Approval Signatures

Written By




R. Kunnen

Reviewed By (optional)




M. Boutelje

Reviewed By (optional)



D. Curiel

Approved By



GJ Plattel

Name: <b>Recycling Passport PFEI-C-30T</b>		12NC: <b>4522 132 20250</b>	Version: <b>04</b> Doc.Date: <b>2007-08-08</b>
Author: <b>R. Kunnen</b>	Pag.: <b>3 of 3</b>	Group No: <b>5790</b>	Section: - Status: <b>Authorized</b>

