




Product name	NT-DAC 3T-C-XW6000 / XW6200 / XW6400
Identification code(s)	4522 132 2013x / 4522 132 5079x / 4522 132 5880x / 9896 030 1859x / 9896 030 1821x
Total weight [kg]	± 137 kg excluding CDAS (±27kg) and XW6000 /6200 /6400 (±17kg)
Manufacturer	Philips Medical Systems Nederland B.V.
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	System_NTDAC3T-XW6000-6200-6400CDAS_452213250790_579_v01.doc

Recycling information

Special Attention	Items	Location
	NOT part of this passport:	
	1/ CDAS (45221328904x / 45221325996x)	Figure 2
	2/ XW6000/6200(45221322734x /45221324545x)	Figure 3
	XW6200/6400 8GB C (45221326520x / ..6519x)	
	3/ Stand. Cab. 33HE + 2HE Lade (4522 131 7443 x)	Figure 1

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Cable bracket front (Zn yellow passivated 5µm)	Figure 2
	(quick acting) Screw (Zn yellow passivated)	Figure 3
	Cable bracket back (Zn yellow passivated 5µm)	Figure 5
	Cable bracket (Zn yellow passivated 5µm)	Figure 6
	Cable relief bracket (Zn yellow passivated 5µm)	Figure 7

Name: Recycling Passport <NT-DAC 3T-C-XW6200>	12NC: 4522 132 50790	Version: 01 Doc.Date: 2007-07-31
Author: <G. Geboers>	Pag.: 1 of 3	Group No: 5790 Section: - Status: Authorized

Locations as mentioned in the Passport (Information in pictures)



Cdas



Stand. Cabinet





Figure 2



Figure 3



XW 6000
XW 6200
XW 6400

Figure 1 Frontview
NT-DAC cabinet




Figure 5




Figure 7



Figure 6



Figure 4 Backview
NT-DAC cabinet



Name: Recycling Passport <NT-DAC 3T-C-XW6200>		12NC: 4522 132 50790		Version: 01 Doc.Date: 2007-07-31	
Author: <G. Geboers>	Pag.: 2 of 3	Group No: 5790	Section: -	Status: Authorized	Template version 2.

Recycling Passport Revision History

Revision	Description	Date
1	Initial version	2007-07-03

Recycling Passport Approval Signatures

Written By G. Geboers	Reviewed By (optional)
Approved By Geert-Jan Plattèl	Approved By

<END OF DOCUMENT>

Name: Recycling Passport <NT-DAC 3T-C-XW6200>		12NC: 4522 132 50790		Version: 01 Doc.Date: 2007-07-31	
Author: <G. Geboers>		Pag.: 3 of 3	Group No: 5790	Section: -	Status: Authorized