




Product name	Philips Micro coils
Identification code(s)	4522 131 81700 Microscopy 23mm 1.5T 4522 131 81710 Microscopy 47mm 1.5T
Total weight [kg]	0.2 Kg
Manufacturer	Philips Medical Systems
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	RFCoils_Philips_Microcoils.doc

Recycling information

Special Attention 	Items	Location
	None	Figure

Batteries 	Items, type	Location
	None	Figure

Hazardous substances 	Items	Location
	Several printed circuit boards inside the coils.	Figures next page
	Beryllium copper contact pins on connectors	Figures next page

Fluids and gases 	Items	Location
	None	

Name Recycling passport Philips micro coils		12NC: 4522 131 81700		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 1 of 3	Group No: 5790	Section: -	Status: Authorized

Locations as mentioned in the Passport (Information in pictures)

Place sketch (es) or photo(s) here, and indicate issues using icons:



Microscopy 23mm



Microscopy 47mm



Connector

Name Recycling passport Philips micro coils		12NC: 4522 131 81700		Version: 00 Doc.Date: 2006-08-04	
Author: Paul Ebben		Page 2 of 3	Group No: 5790	Section: -	Status: Authorized

Recycling Passport Revision History

Revision	Description	Date
00	Initial version	2006-08-04

Recycling Passport Approval Signatures

<p>Written By Paul Ebben</p>	<p>Reviewed By (optional) Geert-Jan Platel Harrie Schellekens</p>
<p>Approved By Paul Ebben</p>	<p>Approved By (optional)</p>

End of document.

Name Recycling passport Philips micro coils	12NC: 4522 131 81700	Version: 00 Doc.Date: 2006-08-04
Author: Paul Ebben	Page 3 of 3	Group No: 5790 Section: -
		Status: Authorized