

Product name	CDAS (PFS-316) (1.5T), CDAS (PFS-316) (3T)		
Identification code(s)	4522 132 0893x 4522 132 0894x	4522 132 5996x 4522 133 0471x	4598 000 1447x 4598 001 8930x
Total weight [kg]	33.6 kg (includes packaging) / 26 kg (without packaging)		
Manufacturer	RadiSys Corporation		
Address	5445 NE Dawson Creek Drive, Hillsboro, OR		
Zip code	97124		
Country	USA		
Electronic info			

Recycling information

Special Attention	Items	Location
	Power Supply, purchased as sealed unit, contents uncertain.	Figure 1


Batteries	Items, type	Location
	CPU board Coin Battery, Lithium	Figure 1

Hazardous substances	Items	Location
	CPU PWB (Remark: PWB = Printed Circuit Board)	Figure 1
	Other PWB's present in other slots (incl. emigaskets)	Not shown
	Fan Controller PWB	Figure 1
	Backplane, Main PWB	Figure 1
	Backplane, R(ear)T(ransission)M(odule) PWB	Figure 1
	EMI gaskets (contact springs) on front- /rearplates / fan tray cover and dustfilter cover contain BeCu	Figure 1
	Power Supply <i>may</i> contain: PWB, Electrolyte Capacitors	Figure 1

Name: Recycling passport CDAS ASSY III	12NC: DMR134332	Version: 00 Doc.Date: 2011-10-05
Author:	Page 1 of 2	Group No: 5790 Section: -
		Status: Approve

Locations as mentioned in the Passport (Information in pictures)

Figure 1: Front of CDAS Chassis Assembly



Name: Recycling passport CDAS ASSY III		12NC: DMR134332	Version: 00 Doc.Date: 2011-10-05	
Author:	Page 2 of 2	Group No: 5790	Section: -	Status: Approve