




<b>Intera</b>	<b>Intera 1.5T Covers</b>
Identification code(s)	CONE FRONT NL ASSY 4522-131-7796x CONE REAR NL ASSY 4522-131-7797x LOWER COVER FRONT 4522-131-7479x LOWER COVER REAR NL 4522-131-7480x SIDE COVER NL (2x) 4522-131-7452x SET TOPCOVERS 4522-131-9903x
Total weight [kg]	± 150 kg
<b>Manufacturer</b>	<b>Philips Medical Systems Nederland B.V.</b>
Address	Veenpluis 4-6 Best
Zip code	5680DA
Country	The Netherlands
Electronic info	SystemCoversIntera_452213177960_579_v01.doc


### Recycling information

Special attention	Items	Location
	Products are made of composite materials consisting of glass, fibers and resins.	
	Microphone printed circuit board in front cover	Figure 1

Batteries	Items, type	Location
	Not applicable	

Hazardous substances	Items	Location
	Not applicable	

Name: <b>CONE FRONT NL ASSY</b>	12NC: <b>4522 131 77960</b>	Version: <b>00</b> Doc.Date: <b>2006-09-11</b>
Author: <b>J. van Vliet</b>	Page <b>1</b> of <b>3</b>	Group No: <b>5790</b> Section: - Status: <b>Auth.</b>

Fluids and gases	Items	Location
	Not applicable	

**Overview of covers:**



Printed circuit board location  
(inside of cover)

**figure 1**

Name: <b>CONE FRONT NL ASSY</b>	12NC: <b>4522 131 77960</b>	Version: <b>00</b> Doc.Date: <b>2006-09-11</b>
Author: <b>J. van Vliet</b>	Page <b>2</b> of <b>3</b>	Group No: <b>5790</b> Section: -
© PHILIPS MEDICAL SYSTEMS NEDERLAND BV, BEST 2006		Status: <b>Auth.</b>

Location

## Recycling Passport Revision History

Revision	Description	Date
1	Initial version	8-9-2006
2	Authorised version	11-9-2006

## Recycling Passport Approval Signatures

<b>Written By</b> J.van Vliet	<b>Reviewed By (optional)</b> G.J. Plattel
<b>Approved By</b> J.van Vliet	<b>Approved By</b> A. Klaver

End of document.

Name: <b>CONE FRONT NL ASSY</b>	12NC: <b>4522 131 77960</b>	Version: <b>00</b> Doc.Date: <b>2006-09-11</b>	
Author: <b>J. van Vliet</b>	Page <b>3 of 3</b>	Group No: <b>5790</b> Section: -	Status: <b>Auth.</b>