





Product name:	7T Achieva MR Scanner Shim Cabinet	
Identification code	4535 673 1972*	
Total weight (in Kg)	Approximately 400 Kg	
Producer /Manufacturer	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	http://www.medical.philips.com/main/company/sustainability/recycling/mri/

Recycle info	Substances:	Location
Hazardous  To be removed	Lead (Pb) (Solder - PC boards)	Figure 1
 Batteries To be removed	Type:	Location
	Not applicable.	
Special attention 	Item	Location
	Storage capacitors.	Figure 1,2
 Fluids / Gases	Item	Location
	Not applicable.	

* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)




Figure 1. 7T Shim cabinet, front view.




Figure 2. 7T Shim Cabinet, rear view.