





## Recycling Passport

<b>Product name:</b>	DigitalDiagnost VM		
<b>Identification code</b>	9890-010-82662		
<b>Order No.</b>			
<b>Device Type</b>	B2B (Business to business)		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
<b>Total weight (in Kg) thereof</b>			<b>421 Kg</b>

Version: VM			
Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	245,303	KG
	iron, high alloy (>5%)	15,245	KG
<b>steel, iron</b>		<b>260,548</b>	<b>KG</b>
nonferrous metals and alloys	aluminum, -alloy	83,993	KG
	copper, -alloy	0,6	KG
	zinc, -alloy	0,092	KG
<b>nonferrous metals and alloys</b>		<b>84,685</b>	<b>KG</b>
glass / ceramics	glass	2,17	KG
	ferrite	0,025	KG
<b>glass / ceramics</b>		<b>2,195</b>	<b>KG</b>
plastics / organic substances	thermoplastic	7,055	KG
	elastomer	5,208	KG
	thermoset	0,582	KG
	lubricating grease	0,022	KG
<b>plastics / organic substances</b>		<b>12,867</b>	<b>KG</b>
batteries	batteries, NiMH	0,036	KG
	batteries, lithium	0,012	KG
<b>batteries</b>		<b>0,048</b>	<b>KG</b>
standard parts	cables	33,679	KG
	power supplies / transformers	5,9	KG
	printed circuit boards	5,529	KG
	brakes / clutches	3,009	KG
	rasters	2,23	KG
	other electronic powered devic	2,105	KG
	wheels / castors	1,98	KG
	mounting parts, attaching part	0,666	KG
	motors, pumps	0,222	KG
	remote controls and similar	0,148	KG
	fans	0,076	KG
<b>standard parts</b>		<b>55,544</b>	<b>KG</b>
relevant materials	lead and -compounds (cables an	5,44	KG
	thallium and -compounds	0,001	KG
	confirmed by supplier		
<b>relevant materials</b>		<b>5,441</b>	<b>KG</b>
prohibited materials	confirmed by supplier		
<b>prohibited materials</b>			
		<b>421,328</b>	<b>KG</b>

Printed copies are uncontrolled documents

## Recycling Passport

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>  <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
	Lead	Fig. 1
	thallium	Fig. 1
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	batteries, lithium	Fig. 2
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
<b>Fluids / Gases</b> 	<b>Item</b>	<b>Location</b>
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport ( Pictures information)

Version: VM



Printed copies are  
uncontrolled documents