





Product name:	<b>7T Achieva MR Scanner RCC Cabinet</b>	
Identification code	<b>9896 030 5053*</b>	
Total weight (in Kg)	Approximately 152 Kg	
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/mri/">http://www.medical.philips.com/main/company/sustainability/recycling/mri/</a>

Recycle info	<b>Substances:</b>	<b>Location</b>
 <p><b>Hazardous</b> <b>To be removed</b></p>	Lead(Pb) (Solder – PC boards)	Figure 1
	Monitor. Possible mercury in display screen. See monitor label.	Figure 1
 <p><b>Batteries</b> <b>To be removed</b></p>	<b>Type:</b>	<b>Location</b>
	2 Lithium button cell batteries.	Figure 1
 <p><b>Special attention</b></p>	<b>Item</b>	<b>Location</b>
	Not applicable	
 <p><b>Fluids / Gases</b></p>	<b>Item</b>	<b>Location</b>
	Not applicable	

\* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)




Figure 1. 7T RCC cabinet, front view without covers.




Figure 2. 7T RCC Cabinet, rear view without covers.