





<b>Product name:</b>	<b>Couch Assembly, 7T</b>	
<b>Identification code</b>	4535 670 5502*	
<b>Total weight (in Kg)</b>	approximately 204 Kg	
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	<a href="http://www.medical.philips.com/main/company/sustainability/recycling/mri/">http://www.medical.philips.com/main/company/sustainability/recycling/mri/</a>

Recycle info	Substances:	Location
<b>Hazardous</b>  <b>To be removed</b>	Lead (Pb) (solder – PC boards)	Figures 1-5
 <b>Batteries</b> <b>To be removed</b>	Type:	Location
	Lead-acid storage battery	Figure 1
<b>Special attention</b> 	Item	Location
	Not applicable	
 <b>Fluids / Gases</b>	Item	Location
	Not applicable	

\* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport ( Pictures information)




Figure 1: Basic Triggering Unit (BTU) beneath rear lower cover




Figure 2: Peripheral Pulse Unit (PPU) inside table handle




Figure 3: FastTrack Patient Communication FTCOM board location under front end



Figure 4: FastTrack Patient Communication (FTCOM) board detail




Figure 5a. Fastrack Controller top view



Figure 5b Fastrack Controller bottom view.