





Product name:	Gemini GXL (PET)	
Identification code(s)	4535 679 75681, 4535 679 75691, 4535 679 71891	
Total weight	2500 lbs	
Producer/ Manufacturer	Name company:	Philips Medical Systems Cleveland
	Address:	595 miner Road, Cleveland OH
	Zip code:	44143
	Country:	United States

Recycle Info	Items:	Location/(Qty)
Special attention 	Gas Springs	Figure 2 (2)
Fluids / Gases 	None	
Batteries  To be Removed	Type: Sealed Lead Acid - PET Gantry UPS (Optional)	Location Figure 10 (1)

Hazardous	Substances:	Location/(Qty)
 <p>To be Removed</p>	Lead (Pb) – Front Inner Ring	Figure 1 (1)
	Lead (Pb) – Front Outer Ring	Figure 1 (1)
	Lead (Pb) – Rear Inner Ring	Figure 1 (1)
	Lead (Pb) – Rear Outer Ring	Figure 1 (1)
	Lead (Pb) – IVS	Figure 1 (1)
	12V Power Supplies	Figure 3 (2)
	Printed Circuit Board Assembly(PCBA) (Pb)	Figure 4 (3)
	Printed Circuit Board Assembly(PCBA) (Pb) – (inside Cardcage)	Figure 5 (12)
	Switching Power Supply	Figure 5 (1)
	Printed Circuit Board Assembly(PCBA) (Pb) – Power Supply	Figure 5 (1)
	Power Supply	Figure 5 (1)
	Printed Circuit Board Assembly(PCBA) (Pb) – Power Supply	Figure 5 (1)
	Computer	Figure 6 (1)
	Power Supply	Figure 7 (1)
	Control Panel	Figure 8 (3)
	Printed Circuit Board Assembly(PCBA) (Pb) – PMTs (inside Detector)	Figure 9 (420)
	Printed Circuit Board Assembly(PCBA) (Pb) – (inside Detector)	Figure 9 (7)
	Printed Circuit Board Assembly(PCBA) (Pb)	Figure 9 (7)
	Transformer	Figure 9 (1)
	Printed Circuit Board Assembly(PCBA) (Pb) – Transformer	Figure 9 (1)
Printed Circuit Board Assembly(PCBA) (Pb) – PET Gantry	Figure 10 (1)	
UPS (Optional)		

Pb
Rear Inner &
Outer Ring



Pb
IVS



Pb
Front Inner &
Outer Ring




Figure 1



Gas
Springs




Figure 2




Figure 3



Pb
PCBA



Power
Supplies




Figure 4



Pb
PCBA



Pb
PCBA
Power Supply




Figure 5




Figure 6



Pb
PCBA
Power Supply



Computer



Figure 7



Power
Supply



Pb
PCBA



Pb
PCBA



Pb
PCBA
Transformer




Figure 8





Figure 9



Pb
PCBA



Sealed
Lead
Acid

PET Gantry UPS
(Optional)

Figure 10