





Recycling Passport

Product name:	Practix Convenio	
Identification code	9890-010-83631	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
Total weight (in Kg) thereof		508,3 Kg

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	cast iron, cast steel	148,5	KG
steel, iron		148,5	KG
nonferrous metals and alloys	aluminum, -alloy	62	KG
nonferrous metals and alloys		62	KG
plastics / organic substances	elastomer	12	KG
	oil	5	KG
plastics / organic substances		17	KG
batteries	batteries, lead	138	KG
	batteries, lithium	0,004	KG
batteries		138,004	KG
standard parts	cables	18	KG
	wheels / castors	14	KG
	motors, pumps	14	KG
	power supplies / transformers	12	KG
	printed circuit boards	9	KG
	other electronic powered devic	8	KG
	X-ray tubes	1,24	KG
	remote controls and similar	0,35	KG
standard parts		76,59	KG
relevant materials	lead and -compounds (cables an	66,2	KG
	confirmed by supplier		
relevant materials		66,2	KG
prohibited materials	confirmed by supplier		
prohibited materials			

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	lead	Fig. 1 (1+2)
Batteries  To be removed	Type:	Location
	Lithium Battery	Fig. 1 (3)
	Lead Battery	Fig. 1 (4)
Special attention 	Item	Location
	Discharge Power Capacitors before handling the unit	Fig. 1 (5)
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1 (1)
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)




Fig. 1

Printed copies are
uncontrolled documents