





Recycling Passport

Product name:	Practix 33 Plus		
Identification code	9890 010 81801		
Order No.			
Device Type	B2B (Business to business)		
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	See Philips Intranet	
Total weight (in Kg) thereof			150 KG

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	cast iron, cast steel	95,6	KG
steel, iron		95,6	KG
nonferrous metals and alloys	aluminum, -alloy	6,5	KG
nonferrous metals and alloys		6,5	KG
glass / ceramics	ferrite	0,1	KG
glass / ceramics		0,1	KG
plastics / organic substances	oil	6,5	KG
	elastomer	2,5	KG
plastics / organic substances		9	KG
batteries	batteries, lithium	0,002	KG
batteries		0,002	KG
standard parts	power supplies / transformers	18,2	KG
	other electronic powered devic	8,05	KG
	wheels / castors	5,5	KG
	printed circuit boards	2,2	KG
	cables	2	KG
	X-ray tubes	0,8	KG
	remote controls and similar	0,3	KG
	compensation filters	0,25	KG
standard parts		37,3	KG
relevant materials	lead and -compounds (cables an	1,5	KG
	confirmed by supplier		
relevant materials		1,5	KG
prohibited materials	confirmed by supplier		
prohibited materials			

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1 (1)
Batteries  To be removed	Type:	Location
	Lithium Battery (inside IC)	Fig. 1 (2)
Special attention 	Item	Location
	Discharge Power Capacitors before handling the unit	Fig. 1 (3)
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1 (1)
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)




Fig. 1

Printed copies are
uncontrolled documents