





Product name:	7T Achieva MRI TR Brain Coil
Identification code	4535 673 4757*
Total weight (in Kg)	Approximately 9.54 Kg

Producer /Manufacturer	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	http://www.medical.philips.com/main/company/sustainability/recycling/mri/

Recycle info	Substances:	Location
Hazardous  To be removed	Lead(Pb) (Solder – PC boards)	Figure 1,2

 Batteries To be removed	Type:	Location
	Not applicable.	

Special attention 	Item	Location
	Not applicable	

 Fluids / Gases	Item	Location
	Not applicable	

* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)



Figure 1. 7T Achieva TR Brain Coil. Circuit boards internal to housing.



Figure 2. Rear view of TR Brain coil. Note that extensive circuitry is contained in rear backplane of coil.

