





Product name:	7T Achieva MRI TR Interface box (includes Classic and Nova)	
Identification code	4535 674 1440* (Nova) 4535 674 0422* (Classic)	
Total weight (in Kg)	Approximately 2.6 Kg	
Producer/ Manufacturer	Name company:	Philips Medical Systems Cleveland
	Address:	595 Miner Rd., Highland Hts, Ohio
	Zip code:	44143
	Country:	USA
	Electronic info:	http://www.medical.philips.com/main/company/sustainability/recycling/mri/

Recycle info	Substances:	Location
Hazardous  To be removed	Lead(Pb) (Solder – PC boards)	Figure 2
 Batteries To be removed	Type:	Location
	Not applicable.	
Special attention 	Item	Location
	Not applicable	
 Fluids / Gases	Item	Location
	Not applicable	

* Denotes revision level wildcard designation. (Numbers range from 1-9).

Locations as mentioned in the Passport (Pictures information)



Figure 1. 7T Achieva TR Interface box, enclosed view.




Figure 2. 7T Achieva TR Interface box, opened view.