

Recycling Passport

Product name:	MammoDiagnost analog	
Identification code	9890-010-83981	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
Total weight (in Kg) thereof		453 Kg





9890-010-83981	MammoDiagnost family
9890-010-83001	MammoDiagnost Screening
9890-010-83011	MammoDiagnost Diagnosis
9890-010-83021	MammoDiagnost Multi Care
9890-010-83791	MammoDiagnost Screening plus

Version: MammoDiagnost family			
Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	265	KG
	iron, high alloy (>5%)	67	KG
	cast iron, cast steel	14	KG
steel, iron		346	KG
nonferrous metals and alloys	aluminum, -alloy	24	KG
	copper, -alloy	12,5	KG
	molybdenum, -alloy	0,45	KG
	tungsten, -alloy	0,05	KG
nonferrous metals and alloys		37	KG
glass / ceramics	glass	5,7	KG
glass / ceramics		5,7	KG
plastics / organic substances	thermoplastic	13	KG
	oil	0,6	KG
	elastomer	0,3	KG
plastics / organic substances		13,9	KG
standard parts	cables	18	KG
	mounting parts, attaching part	15	KG
	printed circuit boards	7,5	KG
standard parts		40,5	KG
relevant materials	lead and -compounds (cables an	5	KG
	tungsten and -compounds	0,05	KG
	beryllium and -compounds	0,001	KG
	confirmed by supplier		
relevant materials		5,051	KG
prohibited materials	lead and -compounds	5	KG
	confirmed by supplier		
prohibited materials		5	KG
		453,151	KG

Printed copies are uncontrolled documents

Recycling Passport

Printed copies are
uncontrolled documents

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1
	tungsten	Fig. 1
	beryllium	Fig. 1
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1
	HV-transformer oil contains no PCBs	

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Printed copies are
uncontrolled documents

