





Recycling Passport

Product name:	MRS	
Identification code	0000-007-04012	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
Total weight (in Kg) thereof		436 Kg

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
steel, iron	iron, low alloy (<5%)	217,5	KG
	iron, high alloy (>5%)	3,4	KG
steel, iron		220,9	KG
nonferrous metals and alloys	aluminum, -alloy	4,193	KG
	copper, -alloy	2,2	KG
	zinc, -alloy	0,73	KG
nonferrous metals and alloys		7,123	KG
glass / ceramics	ferrite	0,5	KG
	glass	0,062	KG
glass / ceramics		0,562	KG
plastics / organic substances	oil	20	KG
	thermoset	3,552	KG
	thermoplastic	2,139	KG
	lubricating grease	0,02	KG
plastics / organic substances		25,711	KG
batteries	batteries, lead	84	KG
batteries		84	KG
standard parts	power supplies / transformers	20,7	KG
	cables	18,892	KG
	printed circuit boards	16,304	KG
	other electronic powered devic	15,144	KG
	patient tables	9	KG
	mounting parts, attaching part	8,8	KG
	wheels / castors	2	KG
standard parts		90,84	KG
relevant materials	lead and -compounds (cables an	6,58	KG
	beryllium and -compounds	0,007	KG
	selenium and -compounds	0,001	KG
	arsenic and -compounds	0,001	KG
	antimony and -compounds	0,001	KG
	confirmed by supplier		
relevant materials		6,59	KG
prohibited materials	PBBs (polybromobiphenyls)	0,002	KG
	PBDEs (polybrominated diphenyl	0,002	KG
	lead and -compounds in outer s	0,002	KG
	cadmium and -compounds	0,001	KG
	confirmed by supplier		
	mercury and -compounds	0	KG
prohibited materials		0,007	KG
		435,733	KG

Printed copies are uncontrolled documents

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead	Fig. 1
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil tank / Remove oil before dismantling	Fig. 1
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)

Generator



1



Stand (tube;Collimator)



Printed copies are
uncontrolled documents