





## Recycling Passport

<b>Product name:</b>	PCR Eleva Userstation	
<b>Identification code</b>	9890-010-83731	
<b>Order No.</b>		
<b>Device Type</b>	B2B (Business to business)	
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	See Philips Intranet
<b>Total weight (in Kg) thereof</b>		<b>24 Kg</b>

Klasse	Merkmalbezeichnung	Gesamtgewicht	Gewichtseinheit
batteries	batteries, lithium	0,003	KG
<b>batteries</b>		<b>0,003</b>	<b>KG</b>
standard parts	computers and accessories	10,25	KG
	monitors	9,9	KG
	cables	2	KG
	other electronic powered devic	1,96	KG
<b>standard parts</b>		<b>24,11</b>	<b>KG</b>
		<b>24,113</b>	<b>KG</b>

Printed copies are uncontrolled documents

## Recycling Passport


<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>  <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	batteries, lithium	Fig. 1
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
<b>Fluids / Gases</b> 	<b>Item</b>	<b>Location</b>
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport ( Pictures information)

## PCR Eleva Userstation



1

Printed copies are  
uncontrolled documents