

Recycling Passport

Product name:		PCR Eleva S Hi-res Reader
Identification code		4512-201-04301
Producer /Manufacturer	Name company:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Electronic info:	See Philips Intranet
Total weight (in kg)		240kg
thereof		
Steel, iron		
	Iron, low alloy (<5%)	191kg
	Iron, high alloy (>5%)	12kg
Non-iron-metals and alloys and compounds		
	Aluminum, -alloy	0kg
	Zinc, -alloy	0kg
Glass / ceramics		
	Glass materials	1.5kg
	Ferrite	0.1kg
Plastics/not-iron-metals/organic substances		
	Thermoplastic	15.3kg
	Elastomeric bond	0.2kg
Standard parts		
	Electronic and electrical components and assemblies	16.4kg
	Cables	3.5kg

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Leads and compounds	①+②+③+④+⑤
Batteries  To be removed	Type:	Location
	Lithium	⑥
Special attention 	Item	Location
	Fragile – Handle with care	
	Photo multiplier	⑦
	Fluorescent light	⑧

Locations as mentioned in the Passport (Pictures information)

