




## Recycling Passport

<b>Product name:</b>		<b>PCR Corado Reader</b>
<b>Identification code</b>		<b>4512-201-02671</b>
<b>Producer /Manufacturer</b>	Name company:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Electronic info:	See Philips Intranet
<b>Total weight (in kg)</b>		<b>282kg</b>
<b>thereof</b>		
<b>Steel, iron</b>		
	<b>Iron, low alloy (&lt;5%)</b>	<b>223.4kg</b>
	<b>Iron, high alloy (&gt;5%)</b>	<b>16.3kg</b>
<b>Non-iron-metals and alloys and compounds</b>		
	<b>Aluminum, -alloy</b>	<b>0kg</b>
	<b>Zinc, -alloy</b>	<b>0kg</b>
<b>Glass / ceramics</b>		
	<b>Glass materials</b>	<b>1.1kg</b>
	<b>Ferrite</b>	<b>0.1kg</b>
<b>Plastics/not-iron-metals/organic substances</b>		
	<b>Thermoplastic</b>	<b>15.5kg</b>
	<b>Elastomeric bond</b>	<b>0.2kg</b>
<b>Standard parts</b>		
	<b>Electronic and electrical components and assemblies</b>	<b>21.8kg</b>
	<b>Cables</b>	<b>3.6kg</b>

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>  <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
	Leads and compounds	①+②+③+④+⑤
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
	Lithium	⑥
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
	Fragile – Handle with care	
	Photo multiplier	⑦
	Fluorescent light	⑧

**Locations as mentioned in the Passport ( Pictures information)**

