

Recycling Passport

Product name:	OPTIMUS 65 R 1T 480V	
Identification code	9890-000-62791	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	http://staging.medical.philips.com/main/company/sustainability/recycling/
Total weight (in Kg)		269
thereof		

steel, iron	iron, low alloy (<5%)	110,981	KG
	iron, high alloy (>5%)	1,706	KG
steel, iron		112,687	KG
nonferrous metals and alloys	copper, -alloy	3,913	KG
	aluminum, -alloy	1,326	KG
	zinc, -alloy	0,438	KG
nonferrous metals and alloys		5,677	KG
glass / ceramics	glass	0,288	KG
glass / ceramics		0,288	KG
plastics / organic substances	oil	44,25	KG
	thermoplastic	5,557	KG
	thermoset	3,947	KG
	elastomer	0,221	KG
plastics / organic substances		53,975	KG
standard parts	power supplies / transformers	30,5	KG
	other electronic powered devices	27,703	KG
	cables	21,942	KG
	printed circuit boards	13,043	KG
	mounting parts, attaching part	2,816	KG
	remote controls and similar	0,45	KG
	fans	0,438	KG
standard parts		96,892	KG
		269,519	KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
	Lithium chrome cell, 3V	Figure 1 (1)
	Lithium chrome cell, 3V	Figure 1 (2)
Special attention 	Item	Location
	Pb is present in the soldering process of PCBs	Figure 1, 2
Fluids / Gases 	Item	Location
	Transformer Oil, Shell Diala	Figure 1 (3)
	HV-transformer oil contains no PCBs	

Recycling Passport




Figure 1



Figure 2