

Recycling Passport

Product name:	AMPLIMAT CHAMBER BUCKY/CHEST	
Identification code	9890-000-01615	
Order No.		
Device Type	B2B (Business to business)	
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH
	Address:	Roentgenstrasse 24
	Zip code:	22335 Hamburg
	Country:	Germany
	Web info:	http://staging.medical.philips.com/main/company/sustainability/recycling/
Total weight (in Kg) thereof		0,9

steel, iron	iron, low alloy (<5%)	0,166	KG
	iron, high alloy (>5%)	0,012	KG
steel, iron		0,178	KG
nonferrous metals and alloys	copper, -alloy	0,02	KG
nonferrous metals and alloys		0,02	KG
plastics / organic substances	thermoplastic	0,522	KG
	elastomer	0,08	KG
plastics / organic substances		0,602	KG
standard parts	other electronic powered devices	0,048	KG
	printed circuit boards	0,032	KG
	cables	0,012	KG
	mounting parts, attaching part	0,004	KG
standard parts		0,096	KG
		0,896	KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
	Pb is present in the soldering process of printed wiring boards.	Figure 1
Fluids / Gases 	Item	Location

Recycling Passport



Figure 1