





## Recycling Passport

<b>Product name:</b>	<b>ROT 500 / ROT-GS 500 / ROT 504 / ROT-GS 504 SRM xxxx</b>		
<b>Identification code</b>	<b>9874-008-63132</b> <b>9890-000-63901</b> <b>9890-000-63661</b> <b>9890-000-63731</b> <b>9874-008-60122</b> <b>9874-008-47122</b> <b>9890-000-03961</b>	<b>9890-000-03842</b>  <b>9890-000-85002</b>  <b>9890-000-63842</b>	
<b>Order No.</b>			
<b>Device Type</b>	<b>B2B (Business to business)</b>		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	<a href="http://staging.medical.philips.com/main/company/sustainability/recycling/">http://staging.medical.philips.com/main/company/sustainability/recycling/</a>	
<b>Total weight (in Kg) thereof</b>			<b>27,0</b>

Printed copies are uncontrolled documents

steel, iron	iron, low alloy (<5%)	3,1	KG
	iron, high alloy (>5%)	0,7	KG
<b>steel, iron</b>		<b>3,8</b>	<b>KG</b>
nonferrous metals and alloys	aluminium and -alloys	6,8	KG
	copper and -alloys	2,2	KG
	molybdenum and -alloys	0,9	KG
<b>nonferrous metals and alloys</b>		<b>9,9</b>	<b>KG</b>
glass / ceramics	ferrite	0,3	KG
	glass	0,1	KG
<b>glass / ceramics</b>		<b>0,4</b>	<b>KG</b>
plastics / organic substances	oil	4,2	KG
	thermoset	2	KG
	elastomer	0,1	KG
<b>plastics / organic substances</b>		<b>6,3</b>	<b>KG</b>
relevant materials	lead and -alloys	5,5	KG
	tungsten and -alloys	0,1	KG
	beryllium and -alloys	0,001	KG
<b>relevant materials</b>		<b>5,601</b>	<b>KG</b>

# Recycling Passport

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b> 	<b>Substances:</b>	<b>Location</b>
	Lead as x-ray shielding inside housing: 5,5 kg	Figure 1
	Beryllium as x-ray window in vacuum envelope: 1 g	Figure 1
<b>To be removed</b>		
<b>Batteries</b> 	<b>Type:</b>	<b>Location</b>
<b>To be removed</b>		
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
<b>Fluids / Gases</b> 	<b>Item</b>	<b>Location</b>
	Oil: 4,2 kg	
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport ( Pictures information)




Figure 1

Printed copies are  
uncontrolled documents