





Recycling Passport

Product name:	ROT 350 / ROT 351 / DA 90 / DA 350			
	RO xxxx	SRO xxxx	DU xxxx	DU xxxx
Identification code	9890-000-85282 9890-000-85302	9874-007-01122 9890-000-63182 9874-004-21122 9890-000-63751 9874-005-08122 9890-000-85831 9890-000-85841 9874-006-11122 9874-094-00362 9890-000-85851 9890-000-03951	9890-000-85641 9890-000-85621 9890-000-85631 9890-000-85611	9890-000-85081 9890-000-85061 9890-000-85051 9890-000-85021
Order No.				
Device Type	B2B (Business to business)			
Producer /Manufacturer	Company name:	Philips Medical Systems DMC GmbH		
	Address:	Roentgenstrasse 24		
	Zip code:	22335 Hamburg		
	Country:	Germany		
	Web info:	http://staging.medical.philips.com/main/company/sustainability/recycling/		
Total weight (in Kg) thereof				22,5

Printed copies are uncontrolled documents

steel, iron	iron, low alloy (<5%)	2,86KG
	iron, high alloy (>5%)	0,2KG
steel, iron		3,06KG
nonferrous metals and alloys	aluminium and -alloys	4,3KG
	copper and -alloys	2,2KG
	molybdenum and -alloys	0,5KG
nonferrous metals and alloys		7 KG
glass / ceramics	glass	0,2KG
		0,2KG
plastics / organic substances	oil	3,5KG
	thermoset	1,9KG
	elastomer	0,1KG
plastics / organic substances		5,5KG
relevant materials	lead and -alloys	5,1KG
	cobalt and -alloys	0,24KG
	tungsten and -alloys	0,08KG
relevant materials		5,42KG

Recycling Passport

Recycle info		
General Hazards (Desinfection)	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
Hazardous Substances  To be removed	Substances:	Location
	Lead as x-ray shielding inside housing: 5,1 kg	Figure 1
Batteries  To be removed	Type:	Location
Special attention 	Item	Location
Fluids / Gases 	Item	Location
	Oil: 3,5 kg	
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

Recycling Passport

Locations as mentioned in the Passport (Pictures information)




Figure 1

Printed copies are
uncontrolled documents