





## Recycling Passport

<b>Product name:</b>	ROT 1001 / ROT-GS 1001 / ROT-GS 1002 / ROT 1003 / ROT-GS 1003 / ROT-GS 1005 MRC xxx xxxx		
<b>Identification code</b>	9890-000-03682 9890-000-03692  9890-000-63041	9890-000-85132 9890-000-85142	9890-000-86061
<b>Order No.</b>			
<b>Device Type</b>	<b>B2B (Business to business)</b>		
<b>Producer /Manufacturer</b>	Company name:	Philips Medical Systems DMC GmbH	
	Address:	Roentgenstrasse 24	
	Zip code:	22335 Hamburg	
	Country:	Germany	
	Web info:	<a href="http://staging.medical.philips.com/main/company/sustainability/recycling/">http://staging.medical.philips.com/main/company/sustainability/recycling/</a>	
<b>Total weight (in Kg) thereof</b>			70,0

steel, iron	iron, low alloy (<5%)	13,4KG
	iron, high alloy (>5%)	6,6KG
<b>steel, iron</b>		<b>20KG</b>
nonferrous metals and alloys	aluminium and -alloys	13,6KG
	molybdenum and -alloys	6,5KG
	copper and -alloys	3,6KG
<b>nonferrous metals and alloys</b>		<b>23,7KG</b>
glass / ceramics	ceramics	1,6KG
<b>glass / ceramics</b>		<b>1,6KG</b>
plastics / organic substances	oil	8,5KG
	thermoset	4,6KG
	elastomer	0,7KG
<b>plastics / organic substances</b>		<b>13,8KG</b>
relevant materials	lead and -alloys	8,6KG
	cobalt and -alloys	1,2KG
	tungsten and -alloys	0,3KG
	beryllium and -alloys	0,01KG
<b>relevant materials</b>		<b>10,11KG</b>

Printed copies are uncontrolled documents

# Recycling Passport

<b>Recycle info</b>		
<b>General Hazards (Desinfection)</b>	For dismantling activities Treatment Facilities must consider the national requirements. For personnel which can come into contact with contaminated material, preventive measures pursuant to national requirements must be taken into account	
<b>Hazardous Substances</b>  <b>To be removed</b>	<b>Substances:</b>	<b>Location</b>
	Lead as x-ray shielding inside housing: 8,6 kg	Figure 1
	Beryllium as x-ray window in vacuum envelope: 10g	Figure 1
<b>Batteries</b>  <b>To be removed</b>	<b>Type:</b>	<b>Location</b>
<b>Special attention</b> 	<b>Item</b>	<b>Location</b>
<b>Fluids / Gases</b> 	<b>Item</b>	<b>Location</b>
	Oil: 8,5 kg	
	HV-transformer oil contains no PCBs	

Printed copies are uncontrolled documents

# Recycling Passport

Locations as mentioned in the Passport ( Pictures information)

Beryllium as x-ray window  
in vacuum envelope

Lead as x-ray shielding  
inside housing




Figure 1

Printed copies are  
uncontrolled documents