





**Philips Medical Systems / Product Recycling Passport**


<b>Product name:</b>	<b>X-ray Tube Housing Assemblies – DA 200 Series -</b>	
<b>Identification code(s)</b>	<b>DA 200</b>	
	<b>Note: Additional suffix letters may be included as part of the product identifier.</b>	
<b>Total weight (in Kg)</b>	<b>84.0kg – 98.0kg</b>	
<b>Producer/ Manufacturer</b>	Name company:	Dunlee, Division of Philips Medical System
	Address:	555 North Commerce Street
	Zip code:	Aurora, Illinois 60504
	Country:	USA

<b>Recycle Info</b>	<b>Items:</b>	<b>Location</b>
<b>Special attention</b> 	Product is heavy in nature.	
	This product contains high vacuumed sub-assembly.	
	Product is intended to return to the manufacturer for proper disposal.	
	Contact your Dunlee/Philips Organization for proper handling and disposal.	
<b>Fluids / Gases</b> 	<b>Items:</b>	<b>Location</b>
	Contained Mineral Insulating Oil (No PCBs)	Entire product
<b>Batteries</b>  <b>To be Removed</b>	<b>Type:</b>	<b>Location</b>
	Not applicable	
<b>Hazardous</b>  <b>To be Removed</b>	<b>Substances:</b>	<b>Location</b>
	Internal Lead (Pb) lining	Housing Material
	Beryllium (Be)	X-ray Tube Window Material

Title:  X-ray Tube Housing Assemblies – DA 200 Series -

Recycling Passport Number:  453575069075 Rev:  A

# Philips Medical Systems / Product Recycling Passport



## Recycling Passport Revision History

Rev	Description	Date
A	First Passport	Aug 2006

Title:  X-ray Tube Housing Assemblies – DA 200 Series - \_\_\_\_\_

Recycling Passport Number:  \_453575069075\_ Rev:  A