





Philips Medical Systems / Product Recycling Passports

Product name:	X-ray Tube Housing Assemblies – DA/MX 100 Series -	
Identification code(s)	DA100, MX100 CT Max	
Total weight (in Kg)	34.5 kg	
Producer/ Manufacturer	Name company:	Dunlee, Division of Philips Medical System
	Address:	555 North Commerce Street
	Zip code:	Aurora, Illinois 60504
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Product is heavy in nature.	
	This product contains high vacuumed sub-assembly.	
	Product is intended to return to the manufacturer for proper disposal.	
	Contact your Dunlee/Philips Organization for proper handling and disposal.	
Fluids / Gases 	Items:	Location
	Contained Mineral Insulating Oil (No PCBs)	Entire product
Batteries  To be Removed	Type:	Location
	Not applicable	
Hazardous  To be Removed	Substances:	Location
	Internal Lead (Pb) lining	Housing Material

Title: X-ray Tube Housing Assemblies – DA/MX100 Series -

Recycling Passport Number: 453575043575_ Rev: A

Philips Medical Systems / Product Recycling Passports



Recycling Passport Revision History

Rev	Description	Date
A	First Passport	Aug 2006

Title: __ **X-ray Tube Housing Assemblies – DA/MX100 Series** - _____

Recycling Passport Number: _453575043575_ Rev: __A__