





Philips Medical Systems / Product Recycling Passport


Product name:	X-ray Tube Housing Assemblies - CTR 1500 Series -	
Identification code(s)	CTR15XX, CTT15XX	
	Note: XX represents the product identifier ranging from 00 thru 99. Additional suffix letters may be included as part of the product identifier.	
Total weight (in Kg)	35.5kg	
Producer/ Manufacturer	Name company:	Dunlee, Division of Philips Medical System
	Address:	555 North Commerce Street
	Zip code:	Aurora, Illinois 60504
	Country:	USA

Recycle Info	Items:	Location
Special attention 	Product is heavy in nature.	
	This product contains high vacuumed sub-assembly.	
	Product is intended to return to the manufacturer for proper disposal.	
	Contact your Dunlee/Philips Organization for proper handling and disposal.	
Fluids / Gases 	Items:	Location
	Contained Mineral Insulating Oil (No PCBs)	Entire product
Batteries 	Type:	Location
	Not applicable	
To be Removed		
Hazardous 	Substances:	Location
	Internal Lead (Pb) lining	Housing Material
To be Removed		

Title: __ X-ray Tube Housing Assemblies - CTR 1500 Series - _____

Recycling Passport Number: _453575066075_ Rev: __A__

Philips Medical Systems / Product Recycling Passport



Recycling Passport Revision History

Rev	Description	Date
A	First Passport	Aug 2006

Title: X-ray Tube Housing Assemblies - CTR 1500 Series -

Recycling Passport Number: 453575066075 Rev: A