





Philips Medical Systems / Product Recycling Passport

Product name:	X-ray Tube Housing Assemblies – B502 Series -	
Identification code(s)	B502, B502H	
Total weight (in Kg)	89kg	
Producer/ Manufacturer	Name company:	Dunlee, Division of Philips Medical System
	Address:	555 North Commerce Street
	Zip code:	Aurora, Illinois 60504
	Country:	USA


Recycle Info	Items:	Location
Special attention 	Product is heavy in nature.	
	This product contains high vacuumed sub-assembly.	
	Product is intended to return to the manufacturer for proper disposal.	
	Contact your Dunlee/Philips Organization for proper handling and disposal.	
Fluids / Gases 	Items:	Location
	Contained Mineral Insulating Oil (No PCBs)	Entire product
Batteries  To be Removed	Type:	Location
	Not applicable	
Hazardous  To be Removed	Substances:	Location
	Internal Lead (Pb) lining	Housing Material
	Beryllium (Be)	X-ray Tube Window Material

Show locations of materials mentioned on the previous sheet. (Picture information)

Title: X-ray Tube Housing Assemblies – B502 Series -

Recycling Passport Number: 453575072375 Rev: A

Philips Medical Systems / Product Recycling Passport



Recycling Passport Revision History

Rev	Description	Date
A	First Passport	Aug 2006

Title: X-ray Tube Housing Assemblies – B502 Series - _____

Recycling Passport Number: 453575072375_ Rev: A_____